

Product Data Sheet

miro-pretape



Binden
Bandages



Intended use:

Thin, elastic and breathable foam bandage for skin protection and as padding material, for use on intact skin

Application by:

Risk class:

- Annex: VIII - Rule: 1

Durability:

60 Months

Product description:

The ultra-thin, open-pored PU foam bandage is ideal as a skin-friendly underlay for tape bandages, sports bandages and medical fixations. Thanks to its high flexibility, it fits every body contour without wrinkles and ensures comfortable wear without restricting mobility. The open-pored structure ensures excellent air and moisture permeability, reliably preventing heat build-up and skin irritation. Thanks to its soft, supple surface, the bandage effectively protects the skin from friction and pressure points – both in sports use and in physiotherapeutic or clinical applications. It can be easily torn lengthwise and crosswise, quickly applied and removed without leaving any residue. This makes the bandage a versatile, comfortable and practical aid for all situations in which the skin must be protected and a secure hold ensured at the same time.

Indication:

- Skin protection under functional tape bandages (sports and kinesiology tapes)
- Underlay for support and compression bandages
- Underlay for bandaging in sports to reduce friction
- Protection of the skin against pressure points during prolonged or tight bandaging
- Fixation and skin protection during the application of cold/warm compresses
- Underlay for orthoses, splints or temporary immobilising bandages
- Use in physiotherapy for padding during movement techniques or the use of aids
- Prevention of skin irritation during frequent bandaging
- Skin-friendly fixation for sensitive skin
- Use as a protective layer during sporting activities to reduce chafing

item-n°	length*	width	GTIN	SU	
72750	27,00m	7,00cm	04044044002008	10 pc(s).	

* Length unstretched

Basic UDI für miro-pretape: 40440440024FS - SRN: DE-MF-000005970

miro

miro Verbandstoffe GmbH
In den Wiesen 3
51674 Wiehl
Germany

Phone +49 2262 7248-0
Fax: +49 2262 7248-48
E-Mail: info@miro-verbandstoffe.de
Webseite: www.miro-verbandstoffe.com

Issue status: 17.02.2026



SCAN ME